



Gilson Quality Control Report Pipette Head Addendum

Pipette Head for PIPETMAX®

Default Calibration Values



MAX8x20 Pipette Head			
Volume Setting	1 μ L Minimum	10 μ L 50%	20 μ L Nominal
Calibration Volume Adjustment	1000	10000	20000
Calibration Offset	990	8140	16140
MAX8x200 Pipette Head			
Volume Setting	20 μ L Minimum	100 μ L 50%	200 μ L Nominal
Calibration Volume Adjustment	20000	100000	200000
Calibration Offset	1660	8040	15960
MAX4x1200 Pipette Head			
Volume Setting	120 μ L Minimum	600 μ L 50%	1200 μ L Nominal
Calibration Volume Adjustment	120000	600000	1200000
Calibration Offset	1982	9669	19264

Refer to the **Add Pipette Head(s)** section in the *PIPETMAX® 268 User Guide* (part number LT255519) or *PIPETMAX® 278 User Guide* (part number LT3212000) for information about scanning or entering pipette head calibration values.